

ENFit® Medication Administration Resource

During the past few months, Aveanna Medical Solutions has distributed information regarding changes to enteral feeding connections. These changes are being made by pump and feeding supply manufacturers to ensure that formula, medication, or fluids do not accidentally get infused into the wrong place. If a patient accidentally receives a tube feeding formula in the wrong line, it can be harmful, or in some cases, even fatal.

As your dedicated enteral provider, our priority is to ensure you have the supplies you need to safely administer tube feeding formula. Since Aveanna Medical Solutions does not hold a pharmacy license and does not dispense medications, we are not able to supply medication administration accessories or provide clinical advice regarding the administration of medications. Please direct **questions regarding how to administer medications using the ENFit® system to your physician or pharmacist for guidance.** The information provided below may help you with these discussions.

The following accessories may be used to administer medications through a feeding tube. Please consult your physician or pharmacist on the proper use of these supplies.

Low-dose ENFit® Syringes
(Small volume syringes that reduce wastage of medications)



Medication Bottle Adaptors
(Attaches to a medication bottle allowing an ENFit® syringe to attach)



Medication Straws
(Attaches to an ENFit® syringe and can be inserted in a medication bottle)



Syringe Caps
(If medication is prepared in advance, this cap prevents loss of medication.)



Syringe adaptors

(Variety of adaptors that can be used to allow the attachment of a legacy syringe to an ENFit® feeding or medication port)



Attaches a legacy syringe to an ENFit® feeding or medication port.

Manufacturer: AMT
Item Code # TRN203



Attaches a legacy slip tip syringe to an ENFit® feeding or medication port

Manufacturer: AMT
Item Code # TRN201



Attaches to an ENFit® syringe to convert it to a legacy slip tip syringe.

Manufacturer: AMT
Item Code # TRN101

Resources where pharmacies may order ENFit® supplies

- <https://gohcl.com/>
- <https://www.medidose.com/>
- <https://mms.mckesson.com/>
- <https://www.medline.com/>
- <https://avanosmedicaldevices.com/>
- <https://www.kentecmedical.com/>
- <https://www.vescomedical.com/>
- <https://www.gbukgroup.com/>

Resources where consumers can order ENFit® supplies:

- <https://www.medline.com/>
- <https://www.udelivermedical.com/>
- <https://www.medidose.com/>
- <https://www.vitalitymedical.com/>
- <https://www.atcmedical.com/>
- <https://www.amazon.com/>
- <https://www.medicaleshop.com/>
- <https://www.saveritemedical.com/>

If you have feedback regarding the ENFit® connector system, you can contact:

- Global Enteral Device Supplier Association (GEDSA) who has taken the lead in coordinating the implementation of ENFit®. info@gedsa.org; or
- American Society of Enteral and Parenteral Nutrition (ASPEN) which is the leading organization of nutrition support experts. <https://www.cognitoforms.com/ASPEN9/HowCanASPENHelpYou>

For more information regarding ENFit®, please reference the manufacturer resources on our website at <https://www.aveannamedicalsolutions.com/enfit.html>

ATENCIÓN: si habla español, tiene a su disposición servicios gratuitos de asistencia lingüística. Llame al 1-888-255-8360 (TTY: 711).

注意：如果您使用繁體中文，您可以免費獲得語言援助服務。請致電 1-888-255-8360 (TTY: 711)。

Aveanna Healthcare, LLC does not discriminate on the basis of race, color, religion, national origin, age, sex, sexual orientation, gender identity or expression, disability, or any other basis prohibited by federal, state, or local law.

© 2022 Aveanna Healthcare, LLC. The Aveanna Heart Logo is a registered trademark of Aveanna Healthcare LLC and its subsidiaries.

Aveanna Medical Solutions, 2460 E Germann Road | Chandler, AZ 85286
p: 866 883 1188 | [aveannamedicalsolutions.com](https://www.aveannamedicalsolutions.com)